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


The setup is inspired by a Fabry–Pérot interferometer. Fiber coupled laser light leads to a sensor head where 4% of the light is reflected back to the detector. The rest passes the cavity, is reflected and also guided to the detector. The graph shows the intensity on a detector of a constantly moved mirror - the interference signal.

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

Tuning the wavelength also causes a similar interference pattern as moving the mirror. The graph shows the intensity on a detector with linear decreasing wavelength over time. The absorption lines of a acetylene cell are used for precise wavelength calibration.

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
The fiber-based set-up is completed by an internal reference cavity. The graph shows the signal of the internal reference with adjusted wavelength change in color and the two cavities in black. Analyzing the interference pattern, a lot of information is revealed about the measured cavity.

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The interference patterns show a Doppler-like effect over the up and down sweep of the laser wavelength. First graph I_1 shows a constant negative velocity of the target, the second one I_2 a constant movement in the other direction.

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

CB 2.3 - Collimated beam, standard sensor head

Interferometer type	Fabry-Pérot
Working range	20 ... 5000 mm
Spot size ($2w_0$)	2.3 mm
Fiber connectors	FC/PC Mating Sleeves
Dimensions	12.7 mm, length 30 mm
Targets	mirror, retroreflector

CB 2.3 APC - Collimated beam without reference reflex for individual usage


Interferometer type	none - only beam shaping
Working range	20 ... 5000 mm
Spot size ($2w_0$)	2.3 mm
Fiber connectors	FC/APC Mating Sleeves
Dimensions	12.7 mm, length 30 mm
Targets	mirror, retroreflector

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
FF 50 - Focused beam with fixed focal length

Interferometer type	Fabry-Pérot
Focal length	50 mm
Spot size ($2w_0$)	0.5 mm
Fiber connectors	FC/PC Mating Sleeves
Dimensions	12.7 mm, length 30 mm
Targets	mirror, high reflective surface

FA 50-1400 - Focused beam with adjustable focal length


Interferometer type	Fabry-Pérot
Focal length	50 ... 1400 mm
Spot size ($2w_0$)	< 1 mm
Fiber connectors	FC/PC Mating Sleeves
Dimensions	12.7 mm, length 30 mm
Targets	mirror, high reflective surface

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
MI SR50:50 - Michelson sensor head with different splitting ratios SR for different reflecting targets

Interferometer type	Michelson
Working range	20 ... 5000 mm
Spot size ($2w_0$)	2.3 mm
Beam splitting ratio	50:50, 80:20, 90:10
Fiber connectors	FC/APC Mating Sleeves
Dimensions	12.7 mm, length 30 mm
Targets	mirror, retroreflector

MI OR SR50:50 - Michelson sensor head with open reference and different splitting ratios SR


Interferometer type	Michelson
Working range	20 ... 5000 mm
Spot size ($2w_0$)	2.3 mm
Reference reflectivity	> 90%
Fiber connectors	FC/APC Mating Sleeves
Dimensions	12.7 mm, length 30 mm
Targets	mirror, high reflective surface

quDIS

Interferometric distance measurement

Key Features

- Signal stability < 0.05 nm
- 20 ... 5000 mm working distance
- 25 kHz bandwidth
- 1 m/s target speed
- 3 sensor axes, multiple devices

quDIS Specifications

Sensor	
Sensor axes	3
Working distance * ¹	20 ... 5000 mm
Resolution	1 pm
Signal stability * ²	< 0.05 nm
Bandwidth	25 kHz
Max. target velocity	1 m/s
Fiber input connectors	FC/APC Mating Sleeves
Interferometer	
Laser source	DFB laser (class 1)
Laser power	< 400 μW
Wavelength (IR)	1 535 nm
Laser linewidth	<5 MHz
Alignment Laser * ³	
Laser source	Fiber-coupled diode
Laser power	< 1 mW
Wavelength	650 nm
Interfaces * ³	
PC interface	USB 3.0
Digital out	AquadB & HSSL
Connector	HDMI
Signal levels AquadB / HSSL	LVTTTL / LVDS

Applications

- Interferometric distance measurement
- Beam interrupt compensation
- Vibration analysis
- Angular measurement
- Environment analysis

Operation	
Operating systems	Windows, Linux
Supplied software	GUI, DLL, LabView, Python, Command line
Alignment support	Numerical, graphical
Hardware	
Dimensions	440 x 350 x 50 mm
Weight	4 kg
Power consumption	< 30 W, at 90 to 264 VAC

AMU - Ambient measurement unit * ³	
Dimensions	38 mm x Ø58 mm
Connector	RJ45
Cable length	4 m
Weight	75 g
Sampling rate	10 Hz
Temperature	±0.1°C (-5 ... 50°C)
Pressure	±1 hPa (300 ... 1100 hPa, 0 ... 65°C)
Relative humidity	± 2 % (10... 90% RH, 5... 55°C)
AMU accuracy	±1 ppm

*1: sensor head dependent *2: RMS@ 100Hz, 2s, 200mm *3: optional available

Disclaimer: The information contained herein is subject to change without notice. quTOOLS shall not be liable for technical or editorial errors or omissions contained herein.